

# 1. Scope & Planning

## STEP 1 - Scope & Planning | Injection Molding - PEEK Lumbar Interbody Fusion Cage

| Field                 | Detail                                                                    |
|-----------------------|---------------------------------------------------------------------------|
| Process Name          | Injection Molding - PEEK Lumbar Interbody Fusion Cage                     |
| Scope Start           | Raw PEEK Pellet Receipt & Material Verification                           |
| Scope End             | Final Dimensional Inspection & Lot Release                                |
| Device Classification | Class II - 510(k) (21 CFR 888.3070)                                       |
| Standard Applied      | AIAG-VDA FMEA 4th Edition 2019                                            |
| Regulatory Framework  | 21 CFR 820 (QSR)   ISO 13485:2016   ASTM F2026                            |
| Device Description    | PEEK Lumbar Interbody Fusion Cage - Spinal Implant                        |
| Team Members          | Mold Engineer, Quality Engineer, Process Engineer, Supplier Quality, R&D; |
| FMEA Facilitator      | [Your Name] - Quality Engineer                                            |
| FMEA Number           | PFMEA-SPINE-001                                                           |
| Revision              | Rev A                                                                     |
| Date                  | 2026-06-22                                                                |
| Next Review Date      | 6 months post-transfer validation                                         |
| Status                | Draft - Internal Review                                                   |

AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485

## 2. Structure Analysis

### STEP 2 - Structure Analysis | System Hierarchy Breakdown

| System Level        | Element Name                    | Function Contribution                                   |
|---------------------|---------------------------------|---------------------------------------------------------|
| Level 1 - System    | Injection Molding Cell          | Thermoplastic implant manufacturing environment         |
| Level 2 - Subsystem | Injection Molding Press         | Melts and injects PEEK at controlled temp/pressure      |
| Level 3 - Component | Mold Tool (Cavity/Core)         | Defines part geometry, surface finish, and tolerances   |
| Level 3 - Component | Barrel / Screw Assembly         | Plasticizes PEEK; controls melt temperature and shear   |
| Level 3 - Component | Hot Runner / Gate System        | Delivers melt uniformly; controls gate vestige height   |
| Level 3 - Component | Material Handling & Drying      | Ensures PEEK moisture <0.02% before processing          |
| Level 3 - Component | Process Controller (HMI)        | Enforces validated parameter recipe (IQ/OQ/PQ)          |
| Level 4 - Interface | CMM / Vision Inspection Station | Dimensional and visual verification post-mold           |
| Level 4 - Interface | Electronic Traveler (DHR)       | Records lot, material CoA, operator, and process params |

AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485

### 3. Function Analysis

#### STEP 3 - Function Analysis | Intended Functions & Requirements

| Process Element            | Intended Function                                                                                   | Requirement Satisfied                                                              |
|----------------------------|-----------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| Mold Tool (Cavity/Core)    | Form cage geometry to OD/ID and footprint tolerances per drawing; achieve Ra <= 1.6 um on endplates | Device drawing dimensional stack; osseointegration surface spec per ASTM F2026     |
| Barrel / Screw Assembly    | Maintain PEEK melt temperature 380-400 C; prevent thermal degradation of polymer                    | PEEK processing spec; material property retention per ISO 10993-1 biocompatibility |
| Hot Runner / Gate System   | Deliver bubble-free melt shot; gate vestige height <= 0.3 mm above parting line                     | Device drawing gate vestige callout; visual inspection acceptance criteria         |
| Material Handling & Drying | Dry PEEK pellets to moisture <= 0.02% at 150 C for >= 4 hours before molding                        | PEEK processing guideline; prevention of hydrolytic degradation and splay defects  |
| Process Controller (HMI)   | Lock and enforce validated parameter recipe; flag deviation >5% from setpoint                       | Validated manufacturing process (IQ/OQ/PQ); 21 CFR 820.70(i) software control      |
| CMM / Vision Inspection    | Verify footprint, height, endplate geometry, and gate vestige against print tolerance               | Inspection plan IP-SPINE-002; acceptance criteria per device drawing Rev D         |
| Electronic Traveler (DHR)  | Gate operation sequence; capture material lot, CoA link, operator, and process params               | 21 CFR 820.184 Device History Record; ISO 13485 7.5.3 traceability                 |

**AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485**

## 4. Failure Analysis

### STEP 4 - Failure Analysis | Failure Modes, Effects & Causes

| Process Step                     | Failure Mode                                                       | Effect on Patient / Downstream Process                                           | Severity (1-10) | Root Cause                                                                           | Current Prevention Control                                            | Current Detection Control                                       |
|----------------------------------|--------------------------------------------------------------------|----------------------------------------------------------------------------------|-----------------|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------|-----------------------------------------------------------------|
| Injection Molding - Parting Line | Flash / fin exceeding 0.1 mm at parting line                       | Sharp edge risk during implant handling; potential particulate in surgical field | 6               | Insufficient clamp force; worn parting line; mold damage                             | Clamp force setpoint in validated recipe; mold PM schedule            | 100% visual inspection per WI-VIS-003; AQL sampling CMM         |
| Injection Molding - Geometry     | Outer footprint or height OOS (dimensional deviation)              | Implant-endplate mismatch; subsidence risk; non-union - patient harm             | 8               | Mold wear / thermal expansion; incorrect recipe offset; shot-to-shot variation       | Recipe lock in HMI; mold qualification per OQ; SPC on critical dims   | CMM inspection per sample plan SP-SPINE-01; first-article 100%  |
| Material Handling                | Wrong PEEK grade loaded (incorrect molecular weight / filler)      | Incorrect mechanical properties; implant fracture or fatigue failure in vivo     | 9               | Manual material selection; similar pellet appearance; label misread                  | Material labeling SOP; operator training                              | Visual label check only - no barcode/CoA lock (DETECTION GAP)   |
| Injection Molding - Fill         | Short shot / incomplete mold fill                                  | Structural porosity in load-bearing cage; fracture risk under spinal loading     | 7               | Insufficient injection pressure; blocked gate; material viscosity variation          | Injection pressure setpoint and transfer position in validated recipe | Visual inspection - porosity may be subsurface and missed (GAP) |
| Gate / Trim                      | Gate vestige height > 0.3 mm above parting surface                 | Mechanical interference with endplate; surgeon re-work; particulate generation   | 7               | Gate freeze inconsistency; dwell time variation; manual trim variability             | Dwell time in recipe; gate vestige spec on drawing                    | First-article caliper check; periodic sampling only             |
| Injection Molding - Surface      | Endplate surface texture OOS (Ra > 1.6 um or peak pattern missing) | Reduced bone in-growth; long-term instability and re-operation risk              | 8               | Worn or damaged mold texture; contamination on cavity surface; improper mold release | Mold surface inspection at PM; release agent SOP                      | Profilometer sample check per lot - not 100% (DETECTION GAP)    |

**AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485**

## 5. Risk Rating

**AIAG-VDA 2019 Action Priority (AP) -- NOT Traditional RPN**

**AIAG-VDA 2019 replaces RPN (SxOxD) with Action Priority (AP: High / Medium / Low). For spinal implants, an undetected material substitution (wrong PEEK grade) reaching a patient represents a catastrophic fracture risk regardless of how rarely it occurs. AP weights Detection gaps heavily - a D=7 with any S>=9 is automatically High priority.**

**AP Lookup: S9-10+D7-10->H | S9-10+D4-6+O>=4->H | S7-8+D7-10+O>=4->H | All others->M or L**

### STEP 5 - Risk Rating | AIAG-VDA Action Priority Table

| Failure Mode                  | S (Severity 1-10) | O (Occurrence 1-10) | D (Detection 1-10) | Action Priority | AP Rationale                                                  | Legacy RPN (Reference Only) |
|-------------------------------|-------------------|---------------------|--------------------|-----------------|---------------------------------------------------------------|-----------------------------|
| Flash / fin at parting line   | 6                 | 4                   | 3                  | L               | Lower severity with adequate controls                         | 72                          |
| Outer footprint or height OOS | 8                 | 4                   | 5                  | M               | S=7-8 + D=4-6 + O>=4: moderate risk profile                   | 160                         |
| Wrong PEEK grade loaded       | 9                 | 3                   | 7                  | H               | S>=9 + D>=7: undetected critical failure - patient risk       | 189                         |
| Short shot / incomplete fill  | 7                 | 3                   | 2                  | L               | S=7-8 + D<=3: detection adequate for severity level           | 42                          |
| Gate vestige height > 0.3 mm  | 7                 | 4                   | 4                  | M               | S=7-8 + D=4-6 + O>=4: moderate risk profile                   | 112                         |
| Endplate surface texture OOS  | 8                 | 5                   | 7                  | H               | S=7-8 + D>=7 + O>=4: high-severity, frequent, poorly detected | 280                         |

**AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485**

## 6. Optimization

### STEP 6 - Optimization / Corrective Actions | Risk Reduction Plan

| Initial AP | Failure Mode                  | Action Description                                                                                                                                                                                                                          | Action Owner     | Target Date | Revised S | Revised O | Revised D | Revised AP |
|------------|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|-------------|-----------|-----------|-----------|------------|
| H          | Wrong PEEK grade loaded       | Implement material barcode scan at press hopper: eDHR locks out start until scanned lot matches approved CoA on traveler. Certificate of Analysis auto-linked to DHR. Validate per 21 CFR 820.70(i); 100% operator training before go-live. | Quality Engineer | 2026-Q3     | 9         | 2         | 2         | M          |
| H          | Endplate surface texture OOS  | Install inline profilometer at ejection station: 100% Ra measurement, SPC alert at Ra > 1.4 um. Mold surface inspection added to PM checklist at every 5,000 shots. Update Control Plan CP-SPINE-001 and WI-MOLD-002 Rev B.                 | Process Engineer | 2026-Q3     | 8         | 3         | 3         | L          |
| M          | Outer footprint or height OOS | Add 100% vision system check for critical footprint and height at ejection. Non-conforming parts automatically rejected to quarantine bin. SPC charts for cavity-specific dimensional trending.                                             | Quality Engineer | 2026-Q4     | 8         | 2         | 2         | L          |
| M          | Gate vestige height > 0.3 mm  | Add gate vestige height to mandatory first-article CMM program (currently omitted). Increase periodic sampling from 1/lot to 5/lot until Cpk >= 1.67 demonstrated over 30 lots.                                                             | Quality Engineer | 2026-Q4     | 7         | 3         | 2         | L          |

# 7. Results Summary

## STEP 7 - Results & Documentation Summary

### pFMEA Metrics Summary (auto-computed)

|                                   |   |
|-----------------------------------|---|
| Total Failure Modes Analyzed      | 6 |
| Initial High (H) Action Priority  | 2 |
| High AP After Actions Implemented | 0 |
| Medium (M) AP - Initial           | 2 |
| Medium (M) AP - After Actions     | 1 |
| Low (L) AP - After Actions        | 5 |
| Total Corrective Actions Assigned | 4 |
| Actions with Owner Assigned       | 4 |
| Actions with Target Date          | 4 |

### Linked Quality System Documents

| Document Type    | Document Number    | Title / Description                                        | Status    |
|------------------|--------------------|------------------------------------------------------------|-----------|
| Control Plan     | CP-SPINE-001 Rev A | Injection Molding - PEEK Spinal Cage Control Plan          | In Review |
| Work Instruction | WI-MOLD-002 Rev B  | Injection Molding Setup & Operation - PEEK Cage            | Released  |
| Inspection Plan  | IP-SPINE-002 Rev A | CMM & Vision Inspection - Lumbar Interbody Cage            | Released  |
| IQ/OQ Protocol   | IQ-MOLD-003        | Injection Press Qualification - PEEK Processing Parameters | Pending   |
| PQ Protocol      | PQ-SPINE-CAGE-001  | Process Validation - PEEK Cage Injection Molding Transfer  | Planned   |
| Risk Management  | RM-SPINE-2024-02   | Device Risk Management File (ISO 14971) - Fusion Cage      | In Review |
| DFMEA            | DFMEA-SPINE-001    | Design FMEA - PEEK Lumbar Interbody Fusion Cage            | Released  |
| SOP              | SOP-QE-022         | pFMEA Creation and Maintenance Procedure                   | Released  |

Archive Location: [PLM System] -> Projects -> SPINE-CAGE-TRANSFER-2024 -> FMEA -> PFMEA-SPINE-001\_RevA.xlsx

AIAG-VDA FMEA 2019 | 21 CFR 820 | ISO 13485